

Air travel information

IATA/UN 38.3




Airlines and other carriers will normally insist on disconnecting / isolating the battery from the Ypush prior to travel. This is for safety reasons.

This is done as follows:

- 1 Raise the right armrest.
- 2 Unplug the battery connector. **61**
- 3 Fit the grey isolating bung to the battery socket. **62**
- 4 Optionally remove the battery from the Ypush as shown earlier in this manual.
- 5 Close the armrest.

The airline or other carrier may want to store the battery either in the hold or the cabin space.

Full information on the battery type and capacity can be found in the Technical Specification section of this manual. This can also be found by scanning the QR code below and on the battery itself. Alternatively visit the following website: www.ypush.com



Packing up the Ypush

For transport or storage



For transportation and or storage purposes it may be necessary to pack up the Ypush in to its travel bag. To do this, perform the following steps:



- 1 Ensure no occupant or animal/s are in the seat and the Ypush is on level ground.
- 2 Remove the ignition key.
- 3 Unzip the roller bag and stand it so that the roller wheels are on the ground.
- 4 Disconnect the battery by withdrawing the battery connection plug. [49](#)
- 5 Fit the grey isolator bung into the battery charge socket. [50](#)
- 6 Remove the battery pack by pulling the battery retainer button and lift the battery clear of the Ypush battery pack compartment. Take care not to damage the battery connection plug. [51](#)

Warning! Store the battery safely - Be sure to observe the battery storage/charging and use temperatures noted in the specification section of this manual.

- 7 Pack the battery in the battery bag to protect it and stow safely. [52](#)
- 8 Grip the foot rest plates and fold them so that they are vertical.
- 9 Fold down the push handles by compressing the safety catch with one hand and folding down with the other. [53](#) [54](#)
- 10 Grip the seat fabric along the front edge and pull up. This will allow the Ypush to start to fold. [55](#)

- 11 Raise both foot rest assemblies until vertical. [56](#)
- 12 Fold the Ypush fully closed by pushing the armrests towards each other. [57](#)
- 13 Fit the frame securing strap in place. [58](#)
- 14 Roll the Ypush into the bag and zip-up fully. [59](#)
- 15 Re-attach the pull strap at the desired height. [60](#)

Packing up the Ypush

For transport or storage



Transporting the Ypush

Attendants: To transport the Ypush in a car or similar, you may need to lift it. The component weights are recorded in the Technical Specification section of this manual.

Do not attempt lifting unless you are confident that you can do so without injury. Ask for assistance if you are not sure.

Note: Make sure that the batteries and Ypush are safely secured before driving the vehicle. Unsecured objects can become lethal in the event of a road-crash/accident.

Assembling the Ypush following transportation

Warning! Never unload your Ypush on sloping, or uneven ground. It could run down the slope and injure someone. Avoid placing it on soft ground as it might get stuck or tip over.



After removing your Ypush and placing it on level ground, assemble it as follows:

- 1 Check parts for any sign of damage before assembling them.
- 2 Ensure the ignition key is in the off position or removed from the Ypush.
- 3 Follow the “Setting up your Ypush” section of this manual.
- 4 Ensure that the brakes are in the safe position. (Make sure the Ypush cannot move when pushed).
- 5 Drive the Ypush without an occupant in-situ to make sure everything is functioning correctly. Refer to the troubleshooting guide and or your dealer / service agent if necessary.